

Chicken IgY

PRODUCT INFORMATION

Product Code: 326058471
Isotype: Chicken IgY fraction
Contents: 1 mg antibody in 1 mL PBS (pH 7.2)
Storage: Below -20°C; avoid repeated freeze/thaw cycles.

SOURCE

Normal chicken IgY is obtained from egg yolk of chickens that have not been immunized.

SPECIFICITY

IgY

PREPARATION

Purified IgY was prepared from egg yolk by dextran sulphate precipitation and sodium sulphate precipitation.

CAUTION

For research use only.
Not for diagnostic and therapeutic use.

REFERENCES

1. Jensenius J, Andersen I, Hau J, Crone M and Koch C. 1981, Eggs: conveniently packaged antibodies. Methods for purification of yolk IgG. Journal of Immunological Methods, 46: 63-68

Manufactured and sold by

SHINO-TEST Corporation

2-29-14 Oonodai, Sagamihara-shi, Kanagawa 229-0011, Japan

TEL: +81-42-753-0354 / FAX: +81-42-786-8553

E-mail: HMGB_info@shino-test.co.jp